**Manual del  
usuario**  
*User's manual*

# TABATABA<sup>®</sup> V3

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*Do not use the TABATABA® V3 in a non-medical environment.*

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*Do not dismantle or interfere with the internal components.*

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*Do not use the TABATABA® V3 in an explosive environment or in the presence of anaesthetic gases.*

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*The TABATABA® V3 is a screening device. Values equal or close to 0 cannot, in any case, guarantee non-poisoning. Only blood testing and clinical examination can confirm a CO intoxication.*

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*The TABATABA® V3 is not a device intended for the protection of goods or persons. It should not be used for domestic (heating, faulty fireplaces) or industrial purposes.*

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*Use exclusively the power supply and accessories provided with the TABATABA® V3 to ensure its performance and safety.*

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*Do not store the TABATABA® V3 near chemical substances that could alter measurements (cleaning products, disinfectants or deodorizers, solvents, adhesives, perfumes, etc.).*

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*The TABATABA® V3 should not be immersed or splashed with liquid to be disinfected.*

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*The TABATABA® V3 is a screening and prevention device in the fight against smoking.*

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*The Ergotrapp® is a single-use anticontamination trapper intended to link the mouth to the device TABATABA® V3.*

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*The Ergotrapp® is designed to contain a volume of air to enable accurate measurement by the sensor. It includes an antibacterial and antiviral filter to reduce the risk of cross-contamination. Furthermore, the non-return valve ensures that the patient does not inhale during the measurement.*

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*The TABATABA® V3 should only be used by health professionals (medical secretary, nurse, doctor, pharmacist, tobaccologist, lung specialist ...) qualified to interpret results and to beware the hygiene and bacterial disinfection rules.*

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*The TABATABA® V3 should not be used in the presence of respiratory illness and/or chest problems.  
In all cases, a medical opinion is required.*

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*After several years of recognized experience with older generations of devices, the TABATABA® V3 device shows its qualities, both in terms of technical performance (ease of use, precision of measurements), and in terms of screening (quality of tests).*

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*The device allows general practitioners and occupational physicians to assess CO poisoning intended for smokers or non-smokers.*

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*The device has no real therapeutic benefits but acts as an incentive and monitoring of smoking cessation and / or assessment of tobacco dependence with a view to adapting the methods of quitting.*

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*There is no limit to the number of exams performed per patient.*

### 1.6. Incidentes o riesgos de incidentes graves | *Serious incidents or risks*

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*In the event of an incident or risk of a serious incident related to the device, health professionals or users can make a declaration to the competent authorities of the Member State. In all cases, the manufacturer must be notified as soon as possible in order to declare and deal with this case of material vigilance.*

## 2. Información técnica | *Technical information*

### 2.1. Material provisto | *Supplied equipment*

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List of items included in the packaging:

- TABATABA® V3 device
- Power supply (see details in §0)
- USB-C cable
- Information sheet
- Calibration certificate
- Transport cover

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<sup>1</sup>Filtration efficiency tests were conducted by an independent laboratory.

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*The Ergotrapp® must be stored in a dry and climate-controlled environment. Extreme variations in temperature and humidity in the storage area can reduce filtration performance.*

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*The TABATABA® V3 complies with the requirements of the EN 60601-1-2:2015 + A1:2021-07 standard regarding the electromagnetic compatibility of medical devices.*

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*The electronic design of the TABATABA® V3 ensures the immunity of the CO sensor to surrounding electromagnetic disturbances.*

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*Therefore, the presence of radiofrequency devices does not affect the CO measurement consistency of the TABATABA® V3.*

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
*Hold down the power button to turn on the device and display the homepage.*

1. *Header: Return to homepage button and battery charge status indicator*
2. *Expired CO measurement*
3. *Ambient CO measurement*
4. *Measurements history*
5. *Parameters and calibration*

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 ***Press the Return to homepage button at any time to stop the ongoing operation and return to the homepage.***

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

*The expired and ambient CO measurement operations are not available during battery charging, the access buttons will appear grayed out on the home page.*

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
*Press the TEST button on the homepage. Insert the Ergotrapp® into the device insertion tube (see picture in 2.2). Start the test when the patient is ready*

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*The patient shall inhale and then hold his breath during the countdown (see § 3.5 to customize the duration) The patient shall then blow slowly until the disappearance of the instruction.*

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

*The end of the test is displayed on the screen, allowing the result display to be postponed. The result is displayed following a color code (see § 3.5 to setup the color levels).*

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*Measure unit can be changed between:  
 PPM (default setting) : CO parts per million by volume of air  
 %HbCO : Percent Carboxyhemoglobin  
 %fHbCO : Percent Fetal carboxyhemoglobin*

PPM	COHb %	fCOHb %		Cigarrillos
0-5 ppm	0,5-0,8	0,8-1,6	No fumador <i>Non-smoker</i>	0 cigarrillos
5-10 ppm	0,8-1,6	1,6-3,1	Fumador pasivo o esporádico <i>Passive or casual smoker</i>	0 cigarrillos
10-15 ppm	1,6-2,4	3,1-4,7	Fumador habitual <i>Frequent smoker</i>	5-10 cigarrillos
15-25 ppm	2,4-4,5	4,7-7,7	Fumador empedernido <i>Heavy smoker</i>	10-20 cigarrillos
25-50 ppm	>4,5	>7,7	Fumador muy empedernido <i>Very heavy smoker</i>	20 cigarrillos

Fuente: Michel-Henri Delcroix, Impact de la mesure du monoxyde de carbone (CO) sur l'arrêt du tabagisme de la femme enceinte : étude T-CAFE ;Impact du CO expiré sur les foetus en salle de naissance ; Grossesse extra-utérine (GEU) et tabac : lien causal, levier pour la prévention (Michel-Henri Delcroix, Impacto de la medición del monóxido de carbono (CO) en la cesación del tabaquismo en la mujer embarazada: estudio T-CAFE; Impacto del CO exhalado sobre los fetos en la sala de partos; Embarazo ectópico (EE) y tabaco: vínculo causal, influencia para la prevención). Medicina humana y patología Universidad de Limoges, 2022. Francés. NNT: 2022LIMO0068

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
*Press the button to save the result according to one of the 3 options below:*

1. *New patient profile creation*
2. *Select and save the result under an already existing patient profile*
3. *Save the result anonymously with an automatic timestamp (ID)*

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*In order to create a new patient profile please type in the following information: Lastname/Name, date of birth, gender*

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*The measurement of ambient CO levels should only be conducted in addition to an analysis of CO in exhaled air for the purpose of raising awareness about passive smoking.*

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*Any value equal to or close to 0 can in no way guarantee non-poisoning. A blood test and a clinical examination must confirm any form of CO intoxication, including passive exposure.*

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Press the  button on the homepage. Select the ambient CO measurement duration among the 3 options:

1. *Real-time: instant start and manual stop*
2. *30min: instant start and automatic stop after 30 minutes*
3. *Custom: start after manually setting the duration (1 to 99 minutes)*

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
*The ambient CO level is displayed continuously during the measurement.*

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*Press the STOP button to interrupt the measurement at any time (also in case of predefined duration, 30min like custom)*

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Once the measurement is completed (either automatically or manually) PPM results are displayed and can be saved.

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In order to save the results, please enter the name of the measurement location.

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
Press the button on the homepage to access the saved patients or ambient data.

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Select the patient or ambient profile to view the saved CO measurements. Press on the button to delete.

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
Press the button on the homepage to view the device serial number, firmware version and adjust the date/time, language, WIFI deactivation, breath holding duration after inhaling and result color settings.

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Adjust the date/time using the arrows. Two formats are available for the date (DDMMYYYY or MMDDYYYY) and time (12h or 24h). Use the arrows to adjust the countdown duration for the holding breath step after inhaling (see §3.2). Likewise for the CO values from which a orange or red color shall be applied to the result.

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
*The date of the last calibration of the cell can be found in the "Cell Calibration" menu under settings.*

*Paragraph 5.2 specifies the rules for periodic maintenance to be followed.*

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

*The sensor calibration must be performed in a well ventilated room complying to the environmental conditions of use specified in section 0.*

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
*Press the  button on the homepage then "Sensor Calibration" to view the dates of the CO sensor change and calibration before starting the calibration steps.*

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
*The calibration starts with a 10 seconds automatic setup of the 0. The TABATABA® V3 must be disconnected from the CO gas bottle and free from any consumable.*

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Connect the CO calibration gas bottle to the TABATABA® V3 by means of the calibration kit accessories (flow regulator, tube, pipette connector), a pipette and a trapper. Enter the gas concentration stated on the bottle calibration certificate rounded to the nearest unit and press the next button.

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


Open the valve of the flow regulator to let the gas flow into the TABATABA® V3. Set the flow rate to 1L/s and press START. CO gas calibration takes 20 seconds.

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
Adjust the measured CO value to the reference concentration of the bottle previously entered.  
 Save to confirm the calibration.

If the difference between the measured CO value and the reference concentration does not meet the sensor accuracy (see 50) the message "**Sensor Error**" will appear. In such case repeat the calibration steps. Shall the error remains, please contact FIM Medical or an authorized distributor.

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Calibration operations are not available during battery charging; the access button will appear grayed out on the settings page.

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
Select the Wifi network named after the TABATABA® V3 serial number:  
Ex : SN=230001 > SSID Wifi=Tabataba-230001

Enter the Wifi security key reminded on the device information sheet, which depends on the TABATABA® V3 serial number:  
Ex : SN=230001 > Mot de passe Wifi=TT230001!

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
Once connected to the Wifi network, open the internet browser and enter the following address in the address bar to access the Webapp interface of the TABATABA® V3:

<http://tabataba.local>

or  
<http://192.168.4.1>

Factory user names/passwords:  
Administrator role: admin/admin  
User role : user/user

Please refer to section 4.6 for more information on user account management (login credentials).

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The dashboard allows to view different information and general details:

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**Device:** serial number, firmware version

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**CO Sensor:** last change date and calibration

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**Memory:** Last saved measurements and memory usage status.

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The memory fill rate is divided into 3 components: the number of patient profiles, the number of patient measurements, and the number of ambient measurements.

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
### Device activity overview

Number of CO measurements chart. The display scale can be adjusted to show data for the past week, month, year, or the past 5 years.

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*Webapp features for patients' data:*

#### **Search of patient's profiles**

*Filter the data using the input fields in the column headers.*

#### **Overview of the last measurements**

*Select a patient's profile to display a chart with the 7 last measurements*

#### **Edition of the patients' details**

*Select one or multiple patients' profile(s)  
Enter the new patient's details  
(lastname/name/date of birth)  
Click on "Apply to selection"*

#### **Data deletion**


*Select the patient's profile(s) to delete  
Click on "Delete"*

#### **CSV data export**


*Select one or multiple patients' profile(s)  
Click on "Export" to download the data as CSV file*

#### **Patient PDF report edition**

*Select a patient's profile  
Click on "PDF Report" to download the patient's data as a pdf file (see below picture)*

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*Filter the ambient data using the input fields in the column headers.*

#### **Data deletion**

*Select the measurements to delete  
Click on "Delete"*

#### **CSV data export**

*Select one or multiple measurements  
Click on "Export" to download the data as CSV file*

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*The device configuration functions (see §3.5) are also available on Webapp, in addition to the following functions:*

### **Editing the device name**

*Using an administrator ID, click on the pencil icon to the right of the device name to modify it.*


*The device name can be used to customize the device designation and appears on the patient PDF report.*

### **Factory reset**

*Using an administrator ID click on "Factory Reset" to erase the device's memory, settings and restore the original user passwords.*

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*The user management interface allows to create, edit or delete Webapp login credentials (see section 4.1).*

*There are two user levels available:*

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### Edit a password

Select the user to modify.

Enter the new password in the right panel.

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*To ensure the protection of data stored on the TABATABA® V3, it is necessary to change the default passwords in accordance with local regulations in force regarding password length and complexity.*

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*The TABATABA® V3 should be cleaned after each use with a soft cloth damped with 70% isopropyl alcohol or with a bactericide-virucide wipe among the following references approved by FIM Medical:*

Bactinyl® toallitas desinfectantes perfumadas  
 Clorox® Healthcare Bleach  
 Sani-Cloth® Bleach/Plus/HB/AF3  
 Super Sani-Cloth®  
 Formula 409®  
 Virex® Plus  
 Mikrozyd® AF wipes  
 Mikrozyd® Universal wipes premium  
 Oxivir Excel® wipes

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*The TABATABA® V3 should not be immersed or splashed with liquid.*

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*The single-use Ergotrapp® must be changed systematically after each patient.*

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**Warning:** Risk of cross contamination if disposable consumables are reused

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**Warning:** Do not sterilize the device or the accessories.

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**Warning:** Used or dirty accessories must be disposed of with biological waste (DASRI).

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*The TABATABA® V3 must be calibrated every year.*

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*The electrochemical cell must be replaced every 5 years as part of maintenance.*

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*Only FIM Medical and its authorized distributors are entitled to carry out maintenance.*

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**Warning:** In the absence of annual calibration measurement drift can lead to biased results.

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**Warning:** The electrochemical cell lifetime is estimated at 5 years.

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*In accordance with the WEEE directive used electronic devices must be treated apart from household waste. The devices must be sent to specific collection sites (waste disposal depot). For more information, you can contact FIM Medical.*

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*Under the terms of the contractual warranty, only repairs are covered. The warranty will only apply if the normal and customary usage conditions of the device have been followed. During the annual maintenance, a series of preventive operations are carried out; however, the maintenance itself cannot guarantee coverage for any malfunctions that may occur after the maintenance has been performed.*

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*FIM Medical assesses that the lifespan of the TABATABA® V3 is 10 years as long as the user respects the mandatory annual maintenance carried out by FIM Medical or one of its authorized partners. FIM Medical cannot be held responsible for any device performance issues in case the user fails to comply with the maintenance obligations and usage conditions.*

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*If the problem persists, or for any other issues, contact FIM Medical or your authorized distributor.*

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***If you would like more  
information, feel free to  
contact us.***

**fim-medical.com**

+33 (0)4 72 34 89 89

contact@fim-medical.com

FIM Medical | 51 rue Antoine Primat | 69100 Villeurbanne | Francia