V03.00.00
Avril 2023

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FIM Medical

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1.Introduction

1.1.List of supplied equipment

The following elements should be found inside the packaging:

- TABATABA® V2 device
- User manual
- Calibration certificate
- 2 x 1.5V AA batteries
- Transport cover

Caution: The TABATABA® V2 should be returned to After Sales Service (for calibration or other operations) in the original packaging, including all padding. Repairs that come under guarantee will not be covered in the absence of the original packaging.

1.2.Device presentation

The TABATABA® V2 is an expired carbon monoxide (CO) analyzer targeted uniquely for smoking cessation. It is a tool used to estimate the CO intoxication of smokers and non-smokers.

Three operating modes:

1. **Expired mode:** measures the level of a person's expired carbon monoxide (in ppm, %COHb and %FCOHb).
2. **Monitoring mode:** measures and records atmospheric CO in a room (used uniquely for nicotine addiction study).
3. **Alarm mode:** visual and sound alarm warning that pre-set threshold has been reached (used uniquely for nicotine addiction study).

1.3.Technical features

TABATABA® V2 features	
CO measurement range	0 - 500 ppm
Accuracy	± 3 ppm maximum
Sensor type	Electrochemical
Sensor lifespan	5 years
Annual drift	10% maximum
Calibration	Annual
Screen	LCD 39 x 15 (mm)
Storage temperature	32 - 122 °F (0 - 50°C)
Temperature for use	32 - 95 °F (0 - 35°C)
Humidity	75% maximum
Altitude of operation	< 6561.7ft (< 2000 metres)
Power supply	3V DC (2 x 1.5V AA)
Current	20mA maximum
Standard references	NF ISO 2859-1, EN 60601-1, EN 60601-1-2, IEC 60601-1-6, EN 62366-1, EN62366-1, EN ISO 10993-1, EN ISO 10993-5, EN ISO 10993-10, NF EN ISO 14971, EN 62304/A1, ISO20417, EN ISO 15223-1, NF EN ISO 13485
Medical class	Im (rule 12)
Software security class	A
GMDN code TABATABA®	35467
Applied part (pipettes + cases)	Type B

Dimensions in packaging	170 x 70 x 30 mm
Total weight of device	170g

First CE marketed : 2007

Pipette features

Storage temperature	32 - 122 °F (0 - 50°C)
Temperature for use	32 - 95 °F (0 - 35°C)
Standard reference	NF ISO 2859-1, EN 62366-1, EN ISO 10993-1, EN ISO 10993-5, EN ISO 10993-10, NF EN ISO 14971, ISO20417, EN ISO 15223-1, NF EN ISO 13485
Medical class	I
GMDN code pipette	44545

Trapper features

Storage temperature	32 - 122 °F (0 - 50°C)
Temperature for use	32 - 95 °F (0 - 35°C)
Standard reference	NF ISO 2859-1, EN 62366-1, EN ISO 10993-1, EN ISO 10993-5, EN ISO 10993-10, NF EN ISO 14971, ISO20417, EN ISO 15223-1, NF EN ISO 13485
Medical class	I
GMDN code trapper	60657

1.4.Symbols

Serial number labels display the following information:

CE marking directive 93/42/CEE or MDR2017/745 + N° of notified body

CE marking MDR 2017/745

Type BF applied part

Device should not be discarded with general household waste. Return it to the manufacturer for disposal.

Consult accompanying documents

Serial number


Manufacturer identification

Do not reuse. Single-use only.

Lot number

Storage temperature: between 0 and 50°C

Expiration date of use

2.Safety

2.1.General

Caution:

- Do not use the TABATABA® V2 in a non-medical environment.
- Do not dismantle or interfere with the internal components.
- Do not use the TABATABA® V2 in an explosive environment or in the presence of anaesthetic gases.
- The TABATABA® V2 is a screening device. Values equal or close to 0 cannot, in any case, guarantee non-intoxication. Only blood testing and clinical examination can determine intoxication.
- The TABATABA® V2 is not a device intended for the protection of goods or persons. It should not be used for domestic (heating, faulty fireplaces) or industrial purposes.
- Use only the accessories supplied by FIM Medical or authorized distributors.
- If the TABATABA® V2 is not to be used for some time, remove the batteries.
- Do not attempt to recharge the batteries.
- Do not expose the device to vibrations or brutal shocks.
- Do not store the TABATABA® V2 in the vicinity of chemical substances (paint, solvents, etc.) that could lead to false measurements.
- Green, orange and red LEDs are not indicators for hardware function alerts. The colour of the LED only indicates different thresholds of intoxication (see § 3.1.1).
- The TABATABA® should not be splashed, wet or immersed in liquid to be cleaned.
- If the batteries show any sign of leakage, it is essential to remove them in the dedicated reprocessing circuits and to clean up any traces of chemical substances.

2.2.Intended use

The TABATABA® V2 is a screening and prevention device in the fight against smoking. Single-use pipettes are intended to link the mouth to the device TABATABA® V2.

Single-use trappers are intended to block a volume of air to allow appropriate measurement by the chemical cell. In addition, the safety valve ensures that the patient does not breathe in during the measurement.

2.3.Operators

The TABATABA® V2 should only be used by health professionals (medical secretary, nurse, doctor, pharmacist, tobaccologist, lung specialist ...).

Health professionals should have no difficulty operating the device.

Results should be interpreted by a doctor trained in nicotine addiction studies.

The operator should be aware of the hygiene rules and bacterial disinfection.

If in doubt, the health professional should refer to the user manual and/or contact FIM Medical or the distributor.

2.4.Patient population

The patient population that can use this device is:

- ➔ People over 14 years of age

2.5.Medical contraindications

The TABATABA® V2 should not be used in the presence of respiratory illness and/or chest problems.

In all cases, a medical opinion is required.

2.6.Clinical benefits and risks associated with the use of the device

After several years of recognized experience with older generations of devices, the TABATABA® V2 device shows its qualities, both in terms of technical performance (ease of use, precision of measurements), and in terms of screening (quality of tests).

The device allows general practitioners and occupational physicians to assess CO poisoning intended for smokers or non-smokers.

The device has no real therapeutic benefits but acts as an incentive and monitoring of smoking cessation and / or assessment of tobacco dependence with a view to adapting the methods of quitting.

There is no limit to the number of exams performed per patient with the TABATABA® V2.

2.7.Serious incidents or risks

In the event of an incident or risk of a serious incident related to the device, health professionals or users can make a declaration to the competent authorities of the Member State. In all cases, the manufacturer must be notified as soon as possible in order to declare and deal with this case of material vigilance.

3.Modes of Operation and Use

3.1.Presentation

The TABATABA® V2 CO detector works with 2 1.5V AA batteries. To install or replace the batteries please refer to §4.1.

3.1.1.Buttons

Lift the green cover to access the mode selector and the selection and validation buttons.

Navigation and selection buttons depending on selected mode

Mode selector:

- Expired mode
- Monitoring mode
- Alarm mode

On/Off button

Unit change or zeroing

3.1.2.Lights

Visual assessment of value:

- Green light: 0-3 ppm
- Green light & orange light: 4-5 ppm
- Orange light: 6-8 ppm
- Orange light & red light: 9-12 ppm
- Red light: Over 12 ppm


3.1.3.Screen

The TABATABA® V2 screen displays text and 4-digit numeric values. Small black triangles indicate mode of use, calculation mode, and the unit of numerical values displayed.

Mode use indicator

Low battery indicator

3.2.Settings

Press "OK/SETUP" to access TABATABA® V2 set-up in Expired mode. Parameter settings:

- Date and time
- Monitoring time
- Two alarm thresholds


Use the "+" and "-" buttons to choose parameters from the scrolled messages: DATE, MONI and ALAR then press "OK/SETUP". Each of these parameters is set using the "+" and "-" buttons and validated with "OK/SETUP".

- Order of time and date setting: day, month, year, hours, minutes
- Monitoring time setting offers choices of 15mins, 1h, 8h or 24h
- Order of alarm settings: low threshold alarm, then high threshold alarm

3.3.Installation of accessories

3.3.1.Single-use trapper

Check the no-return membrane is present inside the trapper.

Warning: If the trapper is not changed and reused between two patients, there is no risk of cross-contamination (bacterial or viral contamination) because the trap is equipped with a non-return valve. For hygienic reasons, the trappers are single-use. Systematically change the soiled trap between each patient.

Warning: Used or dirty accessories must be disposed of with biological waste (DASRI).

Insert the trapper into the intended position in the housing.

3.3.2.Single-use pipette

Warning: The pipettes are single use. Systematically change the pipette between each patient. If the pipette is not changed and re-used between two patients, a risk of cross-contamination (bacterial or viral contamination) is possible during insertion.

Warning: Used or dirty accessories must be disposed of with biological waste (DASRI).

3.4. Modes of operation

3.4.1. Expired mode

Description

Expired mode measures the level of carbon monoxide exhaled by the patient. It also displays an instantaneous measurement of CO in the atmospheric air.

In expired air, the measurement can be displayed:

- in ppm: parts per million
- in % COHb : percentage of CO trapped in the haemoglobin
- in % FCOHb: percentage of CO trapped in the foetal haemoglobin

The relation between ppm and %COHb was shown in a study by M.J. JARVIS, M. BELCHER, C. VESEY, DCS. HUTCHISON, in Thorax, volume 41, n°11, pages 886-887 of November 1986.

The value displayed in expired mode is the peak value, i.e. the highest. Change the measurement unit (ppm, %COHb, %COHbF) by pressing successively on "Unit/0".

The CO trapper captures the last millilitres of expired air, called alveolar air. For an instantaneous measurement in atmospheric air, remove the CO trapper; the electrochemical sensor is in direct contact with the atmospheric air.

Ppm correspondence / number of cigarettes per day		
0-5 ppm	Non-smoker	0 cigarettes
5-10 ppm	Passive or infrequent smoker	0 cigarettes
10-15 ppm	Frequent smoker	10 cigarettes
15-25 ppm	Heavy smoker	20 cigarettes
25-50 ppm	Very heavy smoker	60 cigarettes

Use


1. Place a new single-use CO trapper into the TABATABA® V2. Check the no-return membrane is present inside the trapper.

2. Take a new single-use pipette and check it for visual defects and foreign bodies. Place the pipette inside the trapper and check it fits tightly into the trapper.


3. Press On/Off button


The device displays a series of messages. Follow the instructions:

When the device displays "Insp" ask the patient to breathe in.

When the device displays "Bloc" ask the patient to hold their breath.

When the green LED lights up (device displays "exp") tell the patient to place the pipette in the mouth and breathe out slowly and completely.

4. Wait approximately 15 seconds, while commenting on the rising CO, to make the patient aware.

Between two measurements in Expired mode

1. Discard the used pipette.

2. Remove the used trapper and directly discard it.
3. Switch off the device.
4. Use a dry cloth to clean the opening to remove humidity.

5. Ventilate the sensor and wait one or two minutes.

6. Reinsert a new single-use trapper. Check the trapper is conform (presence of no-return membrane).
7. Take a single-use pipette and check there are no defects or foreign bodies.

Insert the pipette into the trapper and check it fits tightly into the trapper. The device is ready for a new measurement. Remind the patient not to inhale in the pipette.

3.4.2. Mode Monitoring

Description

This mode records the level of atmospheric carbon monoxide over a time defined in the TABATABA® V2 set-up. Monitoring can be run for a recording time of:

- 1/4 hour
- 1 hour
- 8 hours
- 24 hours

In all cases, a measurement is taken every 1 minute. After each recording, the TABATABA® V2 displays the following values:

- TWA (Time Weighted Average)
- STEL (Short Term Exposure Limit)
- Maximum value measured (C)

Use

1. Do not use CO trapper or pipette.
2. Place the TABATABA® V2 cursor on Monitoring mode (M) and switch on via the On/Off button.
3. While pressing on the Unit/0 button, press the On/Off button.
4. The "Init" message appears, and recording begins.
5. Press the On/Off button again and the TABATABA® V2 displays "RUN".

As soon as recording begins, the TABATABA® V2 performs an initial measurement, then switches off to save the batteries. Every minute, the screen displays "RUN" and the measured value, then goes on standby. The TABATABA® V2 can store two monitorings. If a monitoring has already been recorded, the TABATABA® V2 successively displays the corresponding TWA, STEL and C values. While pressing on the Unit/0 button, press the On/Off button to run a new monitoring.

The TABATABA® V2 can store two monitorings, the TWA, STEL and C values displayed corresponding to the values of the last monitoring recorded. If the storage memory is full, the "Full" message is displayed for a few seconds. A long press on the Unit/0 button deletes all monitorings. The "CLR" message appears, indicating that the memory has been deleted.

A monitoring in progress can be stopped by activating the device (On/Off button). The screen displays "RUN".

Validate by pressing the OK/SETUP button.

3.4.3. Alarm Mode

Description

In this mode, a sound and visual alarm is set off when the preset threshold (during set-up of the device) has been reached.

The TABATABA® V2 has two alarm thresholds:

- Low threshold alarm
- High threshold alarm

When an alarm threshold is reached, the device sets off an intermittent sound and visual alarm (red LED). The alarm signal remains active for at least 10 seconds.

These two thresholds are differentiated by a more rapid frequency signal for the high level.

Use


1. Do not use CO trapper or pipette.
2. Place the TABATABA® V2 cursor to Alarm mode.
3. Press the On/Off button. The TABATABA® V2 displays "ALAR". It then goes to low consumption mode (displayer blanks).
4. A measurement is taken every 20 seconds. A short flash is emitted by the green LED if the batteries have sufficient power.
5. Another press on the On/Off button reactivates the device, which displays "Init" then takes an instantaneous measurement and displays it for 20 seconds. The device then returns to low consumption.

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The TABATABA® V2 CO detector operates with two 1.5V AA batteries.

- Use a Phillips screwdriver to unscrew the battery compartment cover support screw
- Insert the batteries following the polarity indicators
- When changing the batteries, remove the batteries and recycle them
- Replace the battery compartment cover
- Screw in the battery compartment cover support screw

4.2.Cleaning

The device housing should be cleaned after each use with a damp cloth and a bactericide-fungicide solution approved by FIM MEDICAL (see §4.3). The use of spray is not recommended as a misdirected jet could damage the sensor.

The trapper is single-use and should be changed systematically after each patient.

The pipette is single-use and should be changed systematically after each patient.

Warning: Risk of cross contamination if disposable consumables are reused

Warning: Do not sterilize the device or the accessories.

Warning: Used or dirty accessories must be disposed of with biological waste (DASRI).

4.3.List of Generic Barcericidal Fungicidal products

Due to a large number of brands and references of disinfectant wipes on the market, FIM Medical has validated references for its products that do not alter the appearance or resistance of the plastic materials in the shells of its equipment.

For the disinfection of its TABATABA® V2, FIM MEDICAL has validated the use of the following wipes or cloths:

- 70% isopropyl alcohol
- Bactinyl® disinfection wipes
- Clorox® Healthcare Bleach
- Sani-Cloth® Bleach
- Sani-Cloth® Plus
- Sani-Cloth® HB
- Super Sani-Cloth®
- Sanicloth® AF3
- Formula 409®
- Virex® Plus
- Mikrozid® AF Wipes
- Mikrozid® Universal Wipes premium

4.4.Maintenance

The TABATABA® V2 must be calibrated annually.

The electrochemical cell must be changed every 5 years during the maintenance.

The calibration expiry date is indicated on the calibration certificate. Please return the device to FIM Medical or an authorized distributor for calibration. Only FIM Medical and its authorized distributors are qualified to perform annual checks of the CO tester. A calibration certificate will be delivered.

Warning: Have the device recalibrated before the expiry date. After this date a drift can falsify results.

Warning: The electrochemical sensor has a lifespan of 5 years.

4.5. Disposal

Conforming to the WEEE directive, electronic devices should not be discarded with general household waste. Devices should be disposed of in specific waste centres. For more information, contact the manufacturer.

4.6. Guarantee

Within the framework of the contractual guarantee, only repairs are covered. The guarantee is only applicable if normal and usual conditions of use are respected. During annual servicing, a certain number of preventive operations are performed; breakdowns following annual service are not covered by the guarantee.

The device is guaranteed for 1 year.

4.7. Service life

FIM Medical estimates that the service life of the TABATABA® V2 is 8 years from the time the user has FIM Medical, or one of its authorised partners, perform the compulsory annual device maintenance.

FIM Medical cannot be held responsible for any faulty device performance if the user does not comply with the maintenance obligations and conditions of use.

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Thank you for reading this manual.
If you would like more information, please
do not hesitate to contact us.